

Exéliea™

Full Life Experience



Hear.

Communicate.

Embrace a Full Life Experience.



PHONAK

life is on



Hear. Communicate. Embrace a Full Life Experience.

Hearing and understanding can sometimes be challenging, especially in difficult and changing listening situations. Thanks to a host of groundbreaking new hassle-free functionalities, Exélia will support your individual hearing needs better than any other available product. With Exélia you can be the very life of the party.



Full performance in a discreet design
that virtually disappears when worn.



Best automatic system

Life is dynamic and no two sound situations are the same. To meet your personal hearing needs, your hearing instruments must constantly adapt to different environments. SoundFlow Premium continually analyzes the sounds around you and automatically, and smoothly, decides which sounds to amplify and which to reduce. This ensures best clarity and comfort in all environments, even when using the phone or watching television.

Optimum hearing performance in loud environments

Probably the biggest challenge is to have a conversation in a noisy environment like a restaurant, in the car or even at home when family and friends are invited. Exélia offers several solutions, one of them being VoiceZoom. VoiceZoom focuses on the speech you want to hear and reduces the annoying background noises. All this is done with a sensational and unmatched frequency resolution of 33 channels!

Clearest sound

For your listening comfort, Exélia is able to detect disturbing facets of the sounds around us like reverberation, wind noise, whistling, sharp impulse sounds, and general environmental noises. For each of these annoying sounds there is a dedicated solution to quickly and expertly eliminate the disturbance. The result: a rich, clear sound experience as natural and as full as life itself.



Be in control of a full life


Exélia is the only instrument that lets you choose the direction you want to hear, a functionality called ZoomControl. Accessed via the myPilot full color display remote control, this feature is paramount in situations like driving a car. Of course, the optional myPilot offers many more possibilities to control your hearing systems.

Full connectivity to the world of modern communication

When combined with the optional iCom universal communication interface, your Exélia hearing instruments easily connect to mobile phones, mp3 players, other Bluetooth enabled devices and FM systems. Listening to telephone conversations, stereo music or television has never been so simple and rewarding. iCom is your personal "all access pass" to the world of communication.



myPilot
command center



iCom
communication interface

Welcome to Exélia – a world without limits

Phonak proudly presents a new paradigm in hearing excellence. Exélia is, without a doubt, the most sophisticated hearing instrument system ever introduced. The unique combination of cutting-edge technology, audiological expertise and wireless connectivity offers you:

- unprecedented hearing performance
- convenient control
- easy access to modern communication and entertainment systems



Exélia celebrates the diversity of personal communication needs with comprehensive, tailored hearing solutions for everyone.

Exélia offers you a Full Life Experience.

