
All the sounds of life

More Hearing. Brilliant Design. CORE Performance.



AUDEO **YES**™

PHONAK

life is on

Audéo YES unveils the full spectrum of sounds – for all your clients

Audéo YES accomplishes what no other hearing system can. The smallest wireless hearing instrument in the world provides the solution to a previously unsolved challenge – Audéo YES broadens the hearing range, restoring high frequencies so the full spectrum of sounds is audible for all your clients. In fact, more hearing is just one of its advantages. Audéo YES is triply unique:

More Hearing

SoundRecover breaks the high frequency barrier. This Phonak exclusive expands audibility, improving speech understanding, as well as awareness of the listening environment, to significantly improve quality of life. By compressing and shifting high frequencies into an area of audible hearing, SoundRecover broadens bandwidth for natural and feedback free reproduction of sound without being harsh or tinny.

Brilliant Design

Audéo YES both exceeds performance demands and perfectly matches individual lifestyles. Its improved, robust Audéo housing is available in 17 fresh color combinations and is consistently chosen by end-users as the most attractive external receiver device. Audéo YES also boasts two external receivers offering the best esthetics and comfort for mild to severe hearing losses.

CORE Performance

Audéo YES, the Canal Receiver Technology (CRT) solution of the CORE Collection, is based on the most sophisticated audio-processing platform in the industry. CORE, the Communication Optimized Real-audio Engine, provides unmatched audiological features, including industry benchmarks for automatic functionality, feedback elimination, directional beamforming and hearing focus control – not to mention full wireless functionality.

Your clients will love the look and benefit from the unparalleled performance, exclusive features and wireless functionality – and so will you. Amplify your portfolio and exceed your clients' expectations with Audéo YES. Whatever their needs and wishes are – the answer is Audéo YES.



The answer is Audéo **YES**



Can I give my clients back the full spectrum of sounds?

YES

SoundRecover is the only effective way to extend hearing by restoring high frequency sounds which would otherwise be inaudible. You will soon consider this exclusive Phonak technology to be an indispensable key to customer satisfaction.



SoundRecover – breaking the high frequency barrier

Audibility across a broad range of frequencies is vital for optimum hearing and understanding. However, receiver and coupling limitations mean that no system can provide sufficient gain of very high frequencies for anything but a mild hearing loss. Furthermore, amplification of such sounds can result in tinny, metallic and harsh sound quality, especially for those trying amplification for the first time.

Another approach to expanding the range of audibility is frequency compression. SoundRecover, the exclusive Phonak non-linear frequency compression system, was first introduced in Naida as a solution for severe to profound hearing losses. By personalizing the amount of compression, it is now possible to apply this proven technique with hearing losses of varying degrees and configurations.

- Optimally suited for sloping hearing losses where high frequency amplification alone is simply not sufficient
- Expands the range of audibility for flatter or gently sloping hearing losses for additional audibility not possible with conventional hearing instruments

From now on, your clients can rediscover all the subtle sounds that complete life's rich soundscapes, naturally and without feedback. Audibility of the full spectrum of sounds improves speech understanding, awareness of the sound environment and listening pleasure – even in difficult listening situations. With SoundRecover, Audéo YES enhances quality of life by reconnecting your clients with the world and truly bringing back all the sounds of life.







Your clients can rediscover all
the subtle sounds that complete
life's rich soundscapes.


Can I meet all my clients' needs with a single device?

YES

Canal Receiver Technology (CRT) offers tremendous flexibility and convenience. In addition to the benefits of smaller size, multiple receivers and coupling options mean you can provide an instant or custom fitting for clients with a wide range of audiometric configurations.


Wide fitting range with two external receivers

Accommodate mild to severe hearing losses simply by choosing the CRT receiver that is most appropriate. Boasting improved positioning and a double cerumen guard, the new receivers offer the best combination of esthetics and comfort.



Standard CRT Receiver

For mild to moderately-severe hearing losses



Power CRT Receiver

For moderate to severe hearing losses

Hassle-free simplicity

From fully open to fully closed acoustic coupling, the complete range of domes and custom shells offers you the chance to provide a comfortable instant or custom fitting solution for every client.



Open Dome



Closed Dome



Power Dome



SlimTip



cShell

Enhanced reliability

All Audéo YES housing components are nano-coated prior to assembly. Consequently, Audéo YES resists sweat and other moisture, guaranteeing stable long-term performance. That is not all: the newly designed battery door ensures a tight seal and improved user ergonomics while light-weight, extra-strong titanium housing pins provide easy use with maximum stability.



Can a hearing system look great with no compromise between size and performance?

YES

With Audéo YES, there is no need to choose between size and functionality. Design and performance merge to offer you a product that satisfies all your clients' hearing needs and cosmetic sensibilities.

Award-winning design – appealing, discreet and intelligent

Audéo YES is housed in the Audéo case, chosen by consumers as the most appealing. Audéo YES is not only beautiful – it also enhances functionality. Audéo YES is incredibly small yet still accommodates a long-lasting 312 battery. Moreover, the optimum microphone alignment guarantees maximum directionality benefits.

- Audéo YES is attractive when shown yet invisible when worn
- The tiny case and lightweight, transparent wire offer unrivalled comfort and discreetness
- Best possible intelligibility in noise is achieved by horizontal microphone-positioning



reddot design award



Back in Black



Lunar Eclipse



Precious Metal



Classic Caddy




Snowblade



Taupe

All the colors of life

Audéo YES is available in 17 color combinations. Color skins can be switched easily to create a unique look. The unusual shape and fun colors also make Audéo YES the perfect choice for older children and teenagers.



Can I guarantee clarity and comfort everywhere?

YES

Audéo YES is based on CORE – the fastest, highest capacity audio-processing platform in the industry. CORE guarantees best performance with the most comprehensive collection of features and options.



SoundFlow – automatically ahead

A noisy restaurant, quiet time with a loved one, a conversation in the chaos of city traffic, her first piano recital. What do these situations have in common? What makes them different? SoundFlow instantly determines this and uses the information to provide the optimum combination of comfort and clarity.

- Allows individual fine-tuning of up to four base programs, setting the boundaries within which SoundFlow functions
- SoundFlow blends these base programs, adaptively activating and integrating the other features at optimum levels

SoundFlow, an automatic system that is two generations ahead of the single-base systems of all competitors, takes automatic functionality to a whole new level. Integrating comfort and clarity enhancing features into a program optimized for each unique listening situation, SoundFlow seamlessly adapts to the ever-changing soundscape.



Can I offer real solutions for difficult listening situations?



VoiceZoom – multi-channel adaptive directional beamforming

VoiceZoom sets the industry benchmark in directional beamforming with the highest frequency resolution in up to 33 channels.

- Automatically and intelligently adapts polar patterns with an optimized directivity index in each beamforming channel
- Allows precise suppression of multiple and moving noise sources

VoiceZoom, the unrivalled solution for improving understanding in background noise, is further enhanced by optimum horizontal microphone alignment achieved by the unique shape of Audéo YES.



EchoBlock System – reverberation cancellation

Indoor reverberation and echoey acoustics diminish clarity by blurring and distorting sounds. Such distortions are no longer a problem with EchoBlock System.

- Instantly removes the disturbing reflections of the reverberant signal
- Sound remains natural and undistorted

EchoBlock System means a better listening experience in rooms with reflective surfaces, large assembly halls, places of worship and other reverberant environments.

NoiseBlock Processing – high resolution noise cancellation

NoiseBlock Processing is particularly effective at reducing uncomfortable environmental noises such as industrial equipment, engines and fans.

- Continuously analyzes the environment, recognizing the subtle differences between speech and non-speech sounds
- Automatically reduces noise, based on the input characteristics of the signal, in up to 20 separate channels, only in those channels where it is needed

Comfortable listening without compromising speech integrity is guaranteed, no matter what the setting.

SoundRelax – impulse noise suppression

Sharp, abrupt sounds like clanging dishes and slamming doors are an obvious irritation, made even worse when amplified. Fast-acting SoundRelax is the effective solution.

- Instantly and accurately identifies impulse sounds
- Applies gain reduction relative to the original peak signal characteristics without altering MPO

SoundRelax cushions otherwise unpleasant impulse sounds without affecting speech clarity or causing unnatural perception of other sounds.



Can my clients enjoy hassle-free listening?

YES

WindBlock Management – wind noise reduction

Wind over a microphone causes unpleasant static and distortion – WindBlock Management effectively reduces this unwanted wind noise.

- Identifies and suppresses noise caused by wind turbulence
- Instruments can remain in directional mode while still effectively managing wind noise

With WindBlock Management, participation in outdoor activities is more pleasant, so your clients can really enjoy the great outdoors.



Real Ear Sound – natural ear acoustics

When microphones are placed outside the pinna, natural front/back and vertical localization cues can be lost. Real Ear Sound solves this problem.

- Applies discrete, frequency-specific directionality to restore the effect of the pinna
- Effectively maintains accurate localization

Correct orientation completes the acoustic sound scene.

ZoomControl – user-controlled hearing focus

ZoomControl is a truly unique feature of the CORE audio-processing platform. No other system in the world can offer clients the ability to choose in which direction to focus hearing.

- Combines real-audio streaming with directional beamforming for choice of hearing focus direction
- Allows focus even to the back or side
- Accessed via myPilot command center

Finally you can offer clients a real solution for situations where they cannot face the speaker, such as when riding a bicycle or conversing in a car.




WhistleBlock Technology – freedom from feedback

In the past, open fittings and power fittings have both been among the configurations most prone to feedback problems. Thanks to the industry benchmark feedback management system, WhistleBlock Technology, those days are over.

- Offers more dynamic stable gain across a wider frequency range than any competitive system
- Differentiates true feedback from naturally occurring tonal signals like those in music so sound quality remains excellent and free of artifact

Besides the obvious benefit of less annoyance for your clients, WhistleBlock Technology also improves fitting flexibility, allowing greater headroom and larger venting for less occlusion.



Now everybody can
have full access to the
wireless world.



Can a small CRT device offer full wireless functionality?


YES

Leading an active life means interfacing with communication and entertainment systems like telephones, televisions, MP3 players, computers or other audio sources. Unlike competitive products, which either have no wireless capability or increase size to achieve it, Audéo YES makes no compromises. It is the smallest external receiver instrument in the world offering full wireless functionality.

CableFree Fitting – iCube offers freedom and flexibility

Harnessing the wireless capabilities of CORE, iCube fulfills a long-standing dream to simplify and enhance the fitting process.

- No more tangled cables, connectors or battery pills
- Fast, smooth and reliable
- Range of up to 10 meters/30 feet between computer and iCube so clients can leave the fitting room and experience real-life sounds while you make adjustments in real time
- Dedicated charger so iCube is easy to store and find and is always ready to use
- Available for Audéo YES and all other members of the CORE Collection in all segments





Wireless freedom

Audéo YES accesses the CORE Hearing instrument Body Area Network (HiBAN) which opens up the world of digital multi-media interaction.

iCom – communication interface

iCom is the communication interface at the center of the HiBAN.

- Wireless connection of Audéo YES to a wide variety of external devices
- Ergonomic design with just one button
- Hands-free, binaural telephoning
- Stereo audio streaming via bluetooth or audiojack
- FM availability with just one FM receiver
- EasyAudio and EasyBluetooth programs automatically activated and separately fine-tuneable for optimized phone usage and audio streaming

myPilot – command center

myPilot offers convenient wireless control, useful status information and more.

- Remote control over program and volume of Audéo YES, including separate left/right adjustments
- ZoomControl access for user-controlled, focused listening in any direction
- Read-out of system status, including battery status, current program and volume levels
- Convenient time/date and alarm functions
- Full-color display
- Customizable program names and icons for easy program identification
- Two attractive color options, high-tech gray or white



Can a hearing system offer CORE benefits for all budgets?



Financial considerations play an important role in product selection. Audéo YES is available at three price levels. While Audéo YES III proves that quality does not have to cost the world, Audéo YES V and IX offer additional benefits and flexibility as mid-range and top-end solutions. The Audéo YES models complete the CORE Collection as the CRT counterparts of Exélia, Versáta and Certéna.

Feature overview

Features/Options	Audéo YES IX	Audéo YES V	Audéo YES III
SoundFlow	Premium	Advanced	Standard
SoundRecover	■	■	■
WhistleBlock Technology	■	■	■
Directionality	VoiceZoom (33 channels)	VoiceZoom (12 channels)	Adaptive digital AudioZoom
NoiseBlock Processing	20 channels	16 channels	6 channels
WindBlock Management	■	■	
Real Ear Sound	■	■	
SoundRelax	■		
EchoBlock System	■		
ZoomControl	■		
DataLogging	■	■	■
Self Learning	■	■	
iCom/myPilot	■	■	■
CableFree Fitting (via iCube)	■	■	■



The answer is Audéo **YES**

Audéo YES combines performance and lifestyle. In other words, it not only looks great – it also works great. With its beautifully small design, industry benchmark features and full wireless capabilities, Audéo YES will convince all your clients, irrespective of their lifestyles or budgets – even those with severe hearing losses. If you are still wondering whether the smallest wireless external receiver device in the world could boost your business, the answer is Audéo YES.

This is how Audéo YES will exceed your clients' expectations:

More Hearing

- SoundRecover finally breaks the high frequency barrier, expanding hearing to include the full spectrum of sounds
- SoundRecover restores audibility of crucial, high-pitched speech sounds and subtle environmental sounds
- With SoundRecover, all sounds are audible without being shrill or tinny
- Speech understanding and awareness of the sound environment are improved, even in difficult listening situations

CORE Performance

- CORE enables full wireless functionality for unlimited connectivity and user control
- CORE also offers unmatched features to guarantee comfort and clarity, including industry benchmarks for:
 - frequency compression
 - automatic functionality
 - feedback elimination
 - directional beamforming
 - hearing focus control

Brilliant Design

- Two external receivers enabling a wide fitting range
- Hassle-free simplicity thanks to comfortable instant and custom fitting solutions
- Optimal microphone positioning for maximum directionality benefits
- Incredibly small size is beautiful when shown but invisible when worn
- Clients' favorite design and winner of the Red Dot Design Award
- Great individual choice thanks to 17 color combinations
- Robust housing guaranteeing enhanced reliability



Give your clients back all the sounds of life.
Offer them Audéo YES.



Life is on

We are sensitive to the needs of everyone who depends on our knowledge, ideas and care. And by creatively challenging the limits of technology, we develop innovations that help people hear, understand and experience more of life's rich soundscapes.

**Interact freely. Communicate with confidence.
Live without limit. Life is on.**